


#196.07 - TOWARDS NEW RELEVANT ALZHEIMER'S DISEASE MODELS FOR TARGET VALIDATION AND DRUG TESTING

Louise MINY¹, Ahmad ALLOUCHE², Sonia KRIDI², Julie COLIN², Nicolas VIOLLE², Jessica RONTARD¹, Benoît Guy Christian MAISONNEUVE¹, Thibault HONEGGER¹
1, NETRI, 69007, Lyon, France
2, ETAP-Lab, 54500, Vandœuvre-lès-Nancy, France

BACKGROUND



Over the past decade, no molecules tested in clinical trials to slow or cure neurodegenerative diseases have been brought to market. We present here an essential first step towards the development of innovative organs-on-chip (OoC) models of Alzheimer's disease, in order to elucidate the underlying mechanisms of the disease and to search for new effective therapies. This innovative high-throughput brain-on-chip platform uses:


- Compartmentalized co-culture of glutamatergic and GABAergic neurons derived from human induced pluripotent stem cells (hiPSCs) with the fluidic isolation of NETRI's DualLink MEA device
- The addition of ETAP-Lab's oligomeric forms of amyloid beta₁₋₄₂ (AβO) in one of the three channels (work in progress).

EXPERIMENTAL DESIGN


COMPARTMENTALIZED MEA-CAPABLE OoC DEVICES

Co-culture of hiPSCs neurons:

- hiPSC-derived glutamatergic neurons (BX-0300) in Channel 1
- hiPSC-derived GABAergic neurons (BX-0400) in Channel 3
- Oligomers and/or compounds applied in Channels 1, 2, or 3
- Response recorded in all channels and microchannels



AβO-INDUCED NEURODEGENERATION



Culture of primary rodent neurons:

- Primary cortical neurons in Channel 1
- At Day 21, AβO challenge for 4 days in Channel 3
- Response analysed in Channel 1 and Channel 3: staining β-tubulin/DAPI and imaging with Operetta CLS High Content Analysis System from REVVITY

ETAP-Lab unique know-how:

- Production of stable and soluble oligomers from their human native proteins Aβ₁₋₄₂
- High reproducible neurotoxicity reversed by reference compounds

RESULTS


RELEVANT HEALTHY ALZHEIMER'S MODEL

Expression of markers by semi-automated quantification using NETRI's proprietary software:

- More than 70% of phenotypic markers (Day 21)
- Less than 10% of pluripotency markers (Day 21)

Increase of firing rate in the neuronal networks-on-chip:

- No cell damages after three recordings (up to Day 21)
- Recording of electrophysiological activity from day 14




RELEVANT INJURY ALZHEIMER'S MODEL

Oligomers for more translational value:

- Oligomer-specific antibody (A11) recognized only AβO (Dot-blot assay)
- Fluorescence intensity increased only upon binding AβO not monomers or fibrils (ANS assay)
- SV-AUC study showed only 8% of monomers in mixture oligomers species, the rest is dimers (60%), trimers (20%) and higher molecular weight species (12%)

Localized axonal degeneration:


- Axonal challenging with AβO induced axonal degeneration in dose-dependent manner
- Soma showed no sign of neurodegeneration after local AβO deposits



CONCLUSION & PERSPECTIVES

By combining NETRI's engineering, biological, and digital expertise with ETAP-Lab's expertise in the modeling of neurological diseases on rodent cells and hiPSCs, and in the development and manufacture of neurotoxic oligomers, we have set up a model of Alzheimer's disease-on-chip.

- Fully differentiation and maturation of human neurons on-chip
- Standard Operating Protocol of human glutamatergic and GABAergic neurons co-culture
- Protocol to induce AβO injury on-chip with rodent cells and hiPSC-derived neurons (in progress)



After developing a model of Alzheimer's disease-on-chip, we will focus on the injury model on hiPSCs and the extraction of digital signatures, using our UpLink™ utility software, to finally add reference compounds and compounds of interest.

Our Brain-on-Chip platform will offer pharmaceutical companies and researchers a new model for preclinical studies, enabling them to reproduce complex neuropathological phenomena finely.



This project is funded by the French Government with "France 2030" for the call for proposals "Innovations in biotherapies & bioproduction"

CONTACT

BOOTH #926

Presented on 2023-11-13
Washington DC, USA

www.netri.com
Phone: +33 4 78 23 08 66
Email: contact@netri.com

www.etap-lab.com
Phone: +33 3 83 44 46 35
Email: contact@etap-lab.com

