

NETRI


# #389.06 - A NEW HUMAN BRAIN ORGANOID-ON-CHIP MODEL MEETING BIOLOGICAL AND INDUSTRIAL REQUIREMENTS OF NEUROLOGICAL PRECLINICAL STUDIES



Héloïse Castiglione<sup>1,2,3</sup>, Camille Baquerre<sup>1</sup>, Lucie Madrange<sup>2,3</sup>, Frank Yates<sup>2,3</sup>, Pierre-Antoine Vigneron<sup>2,3\*</sup>, Jessica Rontard<sup>1\*</sup>, Thibault Honegger<sup>1</sup>

<sup>1</sup> NETRI, 69007 Lyon, France, <sup>2</sup> SupBiotech/CEA-IBFJ-SEPIA, 92260 Fontenay-aux-Roses, France, <sup>3</sup> SupBiotech, Ecole d'ingénieurs en Biotechnologies, 94800 Villejuif, France

## BACKGROUND



Cerebral organoids appear as relevant alternative solutions to model human cerebral cellular organization and development in 3D. Nevertheless, implementing organoids for large-scale drug screening requires reproducibility and scalability. To address this challenge, we propose to generate a Brain Organoid-on-Chip platform with:

- A microfluidic technology for recreating cellular microenvironment
- An on-chip cell culture protocol to improve reproducibility of cerebral organoids
- Compatibility with a full process of cerebral organoid characterization

## EXPERIMENTAL DESIGN

### NETRI'S DUPLEX WELL MICROFLUIDIC DEVICE.


Adapted to 3D cell culture.

Two compartments separated by a porous membrane:

- Open well for 3D culture (5x5x5mm; LxWxH)
- Perfusion channel to enhance nutrient diffusion

Adapted to industrial transfer.

- Standardized fabrication procedures
- Compatible with automation processes and adapted to pharmaceutical readouts (microscopy, MEA, HTS)



## OPTIMIZATION OF CULTURE CONDITIONS ON-CHIP.




Fig. 1. Main steps of guided differentiation protocol into cortical organoids (adapted from Xiang et al., 2019) and three timepoints of EBs introduction into Duplex assayed (hiPSCs: human induced pluripotent stem cells, EB: embryoid body).

- Three timepoints of EBs introduction into Duplex tested (Fig. 1): D+2, D+11, D+18
- Three media renewal conditions tested per timepoint (Fig. 2): 1, 2, 3




Fig. 2. Three conditions of medium renewal in Duplex Well.

Controls: - 24-well plate (P24)  
- 96-well plate (P96)

## RESULTS

### ESTABLISHMENT OF A STANDARD PROTOCOL FOR ORGANOID ON-CHIP CULTURE.

Viability and expected morphology up to D+60

Anti-adherence protocol established to limit organoid adhesion

- Timepoint D+2: partial adherence
- Medium renewal #3: partial adherence

Tendency to obtain more reproducible sizes on-chip compared to 24-well controls




Fig. 3: Morphology of cerebral organoids at D+60 in control conditions and on-chip (Duplex), for different EBs introduction timepoints and methods of medium renewal (5X).




Fig. 4: Cerebral organoids growth from D+2 until D+60 of culture, for control organoids and organoids-on-chip with method #1 of medium renewal (Mean ±SEM, n=3 per group). Statistical differences between the groups were determined by simple linear regression (\*\*\*\*p<0.001); D+2 vs D+11, D+2 vs D+18, D+11 vs D+18: ns, P24 grp1 vs grp3: ns, P24 grp1 vs grp2, P24 grp3 vs grp2: \*\*\*\*).

### CHARACTERIZATION OF BRAIN ORGANOID-ON-CHIP.

Transcript level of expression (RT-qPCR)

Condition (D+11 & renewal #1) on-chip most similar to Controls P24 (Fig. 5)

Expected cytoarchitectures at D+60 with characteristic neurogenic areas (rosettes)

Tend to be more reproducible in terms of structural organization compared to Controls P24 (Fig. 6)




Fig. 5: Fold gene expression values in cerebral organoids on-chip and controls for neural progenitors (SOX2), neurons (TUBB3), anti-apoptotic (BCL2), and pro-apoptotic (BAX) markers.




Fig. 6: Immunofluorescence staining of neural progenitors (SOX2), neurons (TUBB3), astrocytes (GFAP), and mature neurons (MAP2) in cerebral organoids cultured at D+60. (Thunder microscope, Leica, objective 10x).

## CONCLUSION & PERSPECTIVES

By combining NETRI's microfluidic device and cerebral organoids, we have set up a Brain Organoid-on-Chip platform meeting the needs of preclinical applications with organoids exhibiting expected characteristic morphology, cell types, cytoarchitectures, and transcript levels of expression at D+60.



Our cerebral organoids on-chip protocol allows:

- Several possible timepoints of EBs introduction and methods of medium renewal on-chip, depending on the experimental design



Future perspectives include further organoid viability characterization, on-chip analyzes (microscopy, MEA) and complexification of the model by adding cell types for blood-brain barrier (BBB) modelling.

Our Brain Organoid-on-Chip platform paves the way for HTS drug screening and toxicological assessments.



Presented on 2023-11-14  
Washington DC, USA

## CONTACT

www.netri.com  
Phone: +33 4 78 23 08 66  
Email: contact@netri.com

BOOTH  
#926

