



Cerebral organoids

- *In vitro* models of human developing brain
- 3D self-organized cell structures
- Lack of reproducibility & transferability to industrial scale [1]



NETRI microfluidic devices

- Controlled fluid flows
- Cellular microenvironment modeling
- Compatible with industrial readouts & automation processes



Brain-Organoid-on-Chip

- Promising technology to facilitate cerebral organoid scalability for preclinical applications [1]
- Establishment of an on-chip culture protocol for cerebral organoids, with improved reproducibility


MATERIALS AND METHODS

NETRI's DUPLEX WELL MICROFLUIDIC DEVICE

Adapted to 3D cell culture:


Two compartments separated by a porous membrane:

- Open well for 3D culture (5x6x4.6 mm)
- Perfusion channel



Adapted to industrial transfer:

- Standardized fabrication procedures
- Compatible with automation processes
- Adapted to pharmaceutical readouts (microscopy, MEA, HTS)



Membrane anti-adherence method:

- To limit organoid adhesion
- Based on fluid mechanics

OPTIMIZATION OF ON-CHIP CULTURE CONDITIONS

Regionalized differentiation into cortical organoids and assayed conditions:




Fig. 1. Main steps of guided differentiation protocol into cortical organoids (adapted from [2]) and three assayed timepoints of EBs introduction into Duplex Well (hiPSCs: human induced pluripotent stem cells, EB: embryoid body). hiPSCs were reprogrammed from BJ fibroblasts (ATCC, CRL-2522).

- 24-well plate (P24): control organoids cultured in conventional support (+ control)
- 96-well plate (P96): control organoids cultured in 150 L (volume in Duplex) (- control)

• Duplex Well:

→ Three timepoints of EBs introduction (Fig. 1): D+2, D+11, D+18

→ Three conditions of medium renewal (Fig. 2): 1, 2, 3




Fig. 2. Three conditions of medium renewal in Duplex Well.

RESULTS

EXPECTED MORPHOLOGIES UP TO D+60




Fig. 3: Morphology of cerebral organoids at D+60 in control and on-chip conditions, for different EBs introduction timepoints and methods of medium renewal (brightfield, 5X).

- Some organoids adhered on the membrane for Timepoint D+2 and for Condition 3 of medium renewal

TENDENCY TO DISPLAY MORE REPRODUCIBLE SIZES AND GROWTH PROFILES COMPARED TO CONTROLS





Fig. 4: Cerebral organoids growth from D+2 until D+60 of culture, for on-chip and control organoids. On-chip organoids exhibit more reproducible intra- and inter-batch sizes and growth evolutions. (Mean   SEM, n=3 per group). Statistical differences between the groups were determined by simple linear regression (***p<0.001), with Duplex D+2 vs D+11, D+2 vs D+18, D+11 vs D+18: ns; P24 grp1 vs grp3: ns; P24 grp1 vs grp2, P24 grp3 vs grp2: ****).

EXPECTED CYTOARCHITECTURES AT D+60 AND TENDENCY TO EXHIBIT MORE REPRODUCIBLE STRUCTURES COMPARED TO CONTROLS



Structural organizations:

- D+11 and D+18 more optimal than D+2
- 1 and 2 similar
- 3 display most optimal patterns of neurogenic areas (rosettes)

Improved reproducibility for on-chip conditions compared to controls:

- Rosettes with expected morphologies (oval shapes) & well-organized patterns (present along organoid borders)

Fig. 5: Immunofluorescence staining of neural progenitors (SOX2), neurons (TUBB3, MAP2), and astrocytes (GFAP) in cerebral organoids cultured at D+60. (Thunder microscope, Leica, objective 10X).

EXPECTED TRANSCRIPTS LEVELS OF EXPRESSION (RT-qPCR)




Fig. 6: Fold gene expression values in control and on-chip cerebral organoids for neural progenitors (SOX2), neurons (TUBB3), anti-apoptotic (BCL2), and pro-apoptotic (BAX) markers.

CONCLUSION

- On-chip cerebral organoids: expected morphology, cell types, cytoarchitectures & RNA expression levels
- Improved reproducibility: size, growth profile & cytoarchitectures (also demonstrated with another hiPSCs line)
- Several possible timepoints of EBs introduction & methods of medium renewal:

	1	2	3
D+2	✗	✗	✓
D+11	✓	✓	✓
D+18	✓	✓	✓

PERSPECTIVES

- Further organoid viability characterization
- Blood-brain barrier modelling
- Paves the way for drug screening & toxicological assessments

