


DUALINK™ ME A

The DualLink™ MEA is a 3-compartment chip for cell co-culture connected by microchannel technology. Due to their micron scale, only cell extensions can grow within the microchannels, leaving the cell bodies within the compartments themselves. The compartments of the DualLink™ MEA and its microchannel *Tunnels* allow a bilateral connectivity and enhanced fluidic isolation, as well as neuronal functional analysis with integrated microelectrode arrays (MEAs) using Axion BioSystems technology.



TECHNICAL SPECIFICATIONS

Surface Areas	Volumes	Formats	Materials
Channel 1 18.8 mm ² (32.9 mm ² with reservoirs)	Channel 1 3.8 μL (160 μL with reservoirs)	Microfluidic chip 3x2 wells QuarterBento™	Microfluidic chip PolyDiMethylSiloxane biocompatible and low compound absorbing
Channel 2 1.2 mm ² (15.3 mm ² with reservoirs)	Channel 2 0.2 μL (156 μL with reservoirs)	4 chips (52.6x34.6x6.2 mm)	NeoBento™
Channel 3 18.8 mm ² (32.9 mm ² with reservoirs)	Channel 3 3.8 μL (160 μL with reservoirs)	SLAS standard 96-well plate (127.8x85.5x17.1 mm)	Polystyrene (1.4 mm thick + refractive index: 1.59)
			MEA Surface PET (125 μm thick refractive index: 1.64) SU8 (5 μm coating) PEDOT-coated gold electrodes (see NeoBento MEA User Guide for details)

Capabilities


- Co-culture & compartmentalization
- hiPSC derived cell
- Axonal transport
- Functional analysis

Applications

- Drug screening
- Innervated skin
- Toxicology
- Virology
- Neuroinflammation
- Neuromuscular junction
- Motor neuron diseases
- Study of the functional activity of neurons

Readouts

- Electrophysiology
- Immunofluorescence
- Live Dead Assays
- Live Staining
- Liquid chromatography
- Mass Spectroscopy
- Lysis cell/supernatant analysis
- ELISA
- Calcium Imaging



NEOBENTO™ MEA

- Abiding by the 96-well microplate standard format, NeoBento™ MEA is a versatile container that enables electrophysiological recording of microfluidic experimentation while being compatible with standard equipment for liquid handling and imaging.



Confidential document

LABYRINTH

By filling the labyrinth with sterile water, the NeoBento™ MEA will securely hold the water needed to maintain, for several days, a humid environment for your cell cultures.

⚠ Empty the labyrinth before placing in Axion's Maestro

TENSION STRUT

It protects the NeoBento™ MEA against the clamping forces of the Axion Biosystems Maestro. Do not remove it

AXION MEA SURFACE

With up to 48 Electrodes per chip, NeoBento™ MEA has integrated Axion Biosystems MEA technology and is compatible with the Axion Biosystems Maestro Pro for 16-chip recordings

BOTTOM COVER

It protects your cultures when manipulating your NeoBento™ MEA, and it is compatible with low magnification imaging (x10, x20).

⚠ Make sure you remove the Bottom Cover before placing inside Axion's Maestro

LID


The lid protects your cultures and maintains a humid environment.

QUARTERBENTO™

The chips of a NeoBento™ are arranged in 4 PDMS QuarterBento™ adhered to the MEA Surface. They each provide multiple microfluidic architectures while optimizing the oxygenation of cells.

MAESTRO EDGE VERSION


The NeoBento™ MEA version for Maestro Edge has Plugs that take up the space where two QuarterBentos™ cannot be recorded.



Plugs allow you, thanks to their own labyrinths, to bring humidity close to the wells on the edges of the remaining QuarterBentos™.

NEOBENTO™ MEA EDGE vs. PRO

Maestro EDGE Version
8 Chips per plate (the top half)



Maestro PRO Version
16 Chips per plate

